



**University of
Nottingham**

UK | CHINA | MALAYSIA



**University of
Nottingham**

Nanoscale and Microscale Research Centre

nmRC CASE STUDY

**SPRAY DRIED DRUG LOADED
POLYMERIC
MICROPARTICLES**

nmRC_CS_01



University of
Nottingham

UK | CHINA | MALAYSIA



University of
Nottingham
Nanoscale and Microscale Research Centre

Spray Dried Drug Loaded Polymeric Microparticles

An Atomic Force Microscopy (AFM) Case Study



- Spray-dried polymeric microspheres have been investigated for use as a controlled release vehicle for active pharmaceutical ingredient (API) administration by intramuscular injection.
- Combinations of biocompatible polymers can be deployed to balance a targeted release profile with appropriate API solubility and stability in-formulation.
- One such example is a API/PLGA/PVP system (Meeus et al, 2012).
- **The surface characterisation of such systems is crucial in understanding how composition and environment affect physiochemical character and performance.**


X-Ray Photoelectron Spectroscopy (XPS)

- Elemental Composition
- Chemical State Information

Sample ID	Theoretical Atomic Composition (%)			Experimental Atomic Composition (%)		
	C 1s	O 1s	N 1s	C 1s	O 1s	N 1s
Pure						
PLGA	57.9	42.1	—	64.8 ± 0.4	35.2 ± 0.4	—
PVP	75.0	12.5	12.5	78.4 ± 0.5	10.5 ± 0.2	11.1 ± 0.5
PLGA–PVP (wt %)						
20/80	71.6	18.3	10.0	68.3 ± 0.4	27.2 ± 0.3	2.9 ± 0.3
40/60	68.2	24.2	7.6	65.0 ± 0.5	34.2 ± 0.6	0.8 ± 0.2
80/20	64.8	30.1	5.1	64.9 ± 0.8	35.1 ± 0.8	—


Time-of-Flight Secondary Ion Mass Spectrometry (ToF-SIMS)

- Depth Profiling (layer-by-layer chemistry).



Scanning Electron Microscopy (SEM)

- High magnification and resolution morphological and topographical characterisation






- X-Ray Photoelectron Spectroscopy (XPS) can quantify the elemental composition of the particle surfaces and provide information on the chemical state.
- Time-of-Flight-Secondary Ion Mass Spectrometry (ToF-SIMS) can qualify the chemical composition and characterise fine changes on a nanometer resolution with depth.
 - **BUT these techniques cannot resolve the surface of an individual microparticle (size range <10 μ m)**
- Imaging with Scanning Electron Microscopy (SEM) provides visualisation on this scale.
 - **BUT offers no distinction of chemical constituents**
- **None of these techniques make no mechanical assessment of the microspheres**




- **Atomic Force Microscopy (AFM)** is a high resolution scanning probe technique that allows visualisation of material surfaces and their topographical features with nanometre resolution, in addition to enabling surface physicochemical, mechanical, thermal or electrical properties to be measured.

SEM Image



AFM 3D Profile




Surface Characterisation down to the
Nanoscale



- AFM can identify the morphological differences between sample surfaces such as microsphere formulation variants.

API/PLGA/PVP Microspheres




Height Sensor 200.0 nm

10/75/15 (%wt)

Homogeneity

Uniformity




Height Sensor 400.0 nm

20/50/30 (%wt)

Morphology

Form, shape, structure



Height Sensor 800.0 nm

30/25/45 (%wt)


Topography

Distribution of features

Observe changes between formulations




- AFM offers the capacity to generate simultaneous 'multi-channel' data.
- Sample topography can be complemented by simultaneous mapping and quantification of mechanical data including adhesion, stiffness etc.



- **Topographical variability distinguished from material differences i.e. API from polymer domains.**

- AFM allows observation of physical and chemical responses at the surface to formulation parameters, and real-time changes to the environment etc.



E.g. AFM can characterise the surface rearrangement of constituent polymers due to a humidity increase.

Rearrangement and phase separation seen at the sample surface.



- AFM is an invaluable tool in the characterisation of the surfaces of micro- and nano-scale materials such as single spray dried polymeric microparticle formulations
- Nanomechanical data identified measurable phase differences representative of distinct API and polymer domains not resolved by other characterisation methods.
- Real-time, '*in-situ*' changes to single particle surfaces in response to environment observed.
- AFM offers a unique capacity to complement broader chemical and visual analysis with nanometre resolution topographical, and simultaneous mechanical characterisation.
- Insight from AFM develops understanding of formulation dynamics and responses that can be linked to performance and used in development



For more details on the work showcased in this case study see the following publications:

Nanoscale Surface Characterization and Miscibility Study of a Spray-Dried Injectable Polymeric Matrix Consisting of Poly(lactic-co-glycolic acid) and Polyvinylpyrrolidone. Joke Meeus, Xinyong Chen, David J. Scurr, Valeria Ciarnelli, Katie Amssoms, Clive J. Roberts, Martyn C. Davies, Guy Van Den Mooter. *Journal of Pharm. Sci.*, 2012, 101, (9) 3473-3485

Surface Characteristics of Spray-Dried Microspheres Consisting of PLGA and PVP: Relating the Influence of Heat and Humidity to the Thermal Characteristics of These Polymers. Joke Meeus, David J. Scurr, Katie Amssoms, Martyn C. Davies, Clive J. Roberts and Guy Van den Mooter. *Mol. Pharmaceutics* 2013, 10, 3213–3224



- We hope the information provided in this case study is of interest.
- If you wish to get in touch with us to discuss any of the information provided, raise a query/concern or provide feedback then please use any of the methods listed below:

nmRC Commercial Services
Nanoscale & Microscale Research Centre
University Park
Nottingham
NG7 2RD

Telephone: +44(0)115 951 5046
Email: nmcs@nottingham.ac.uk
Fax: +44 (0)115 846 7969
Website: www.nottingham.ac.uk/nmrc-commercial