

University of Nottingham

University of Nottingham Nanoscale and Microscale Research Centre

## nmRC CASE STUDY

VISUALISING PROTEIN AND SURFACTANT DISTRIBUTION IN PLGA MICROSPHERES

nmRC\_CS\_02





# Visualizing Protein and Surfactant Distribution in PLGA Microspheres

A Materials Characterisation Case Study



- Polymers are often used in order to maintain controlled release of drug.
- Water-in-oil-in-water (w/o/w) fabricated PLGA microspheres containing lysozyme stabilised by PVA surfactant were analysed using a suite of materials characterisation techniques to investigate a 2 stage release profile observed from a similar system.





## SURFACE CHARACTERISATION

University of Nottingham Nanoscale and Microscale Research Centre

What Influences the release profile?

What is the morphology of the particles?

How consistent is the size and structure?



Where is the protein component?

How is the surfactant distributed?



000 20mm

## SURFACE CHARACTERISATION METHODS

University of Nottingham Nanoscale and Microscale Research Centre







#### **Pre-Sputtering**



#### **Post-Sputtering**



- AFM imaging shows a surface layer of material measuring approximately
  5nm in thickness.
- After sputtering in the ToF-SIMS this layer is gone, confirming it to be the surfactant.



 Multivariate data analysis is a technique useful in reducing the dimensionality of hyperspectral data to simplify its analysis.

 Multivariate Curve Resolution (MCR), a type of multivariate data analysis, was applied to the microsphere ToF-SIMS chemical data to help identify the pure components in the multicomponent system.

 This analysis was undertaken in addition to testing reference samples of the individual components.



#### **Component Distributions (by MCR)**





## **ToF-SIMS SURFACE CHEMICAL IMAGING**



Green – PLGA Blue – PVA Red – Lysozyme



- ToF-SIMS mapping shows a discontinuous PVA surfactant layer.
- Insufficient stabilisation of the microspheres leads to presence of surface pores.
- Majority of lysozyme observed surrounding these pores Leading to a **burst** release.
- Secondarily localised to pore walls Leading to **slower** release.



University of Nottingham Nanoscale and Microscale Research Centre

#### **3D Localisation of Components**

Green – PLGA Red – Lysozyme

> Video Click to play

Note: Video disabled in pdf version





- The bulk analysis of the microspheres by confocal Raman spectroscopy clearly outlines the composition of the microspheres
- The combination of characterisation techniques employed on these microspheres allows for a schematic for the distribution of components to be made:





- Combinations of materials characterisation techniques can generate a wealth of complementary information critical to the understanding sample behaviours that is not available to one technique alone.
- ToF-SIMS and Raman spectroscopy were successfully used to describe the surface and 3D chemical distribution of the surfactant and protein components of a micro particulate formulation.
- AFM and SEM were used in conjunction with ToF-SIMS to confirm the existance of a 4.5 nm thick overlayer of surfactant (not just a surface feature of PLGA).
- In combination this information together helped explain the development of a 2 phase release profile.



University of Nottingham Nanoscale and Microscale Research Centre

#### For more details on the work showcased in this case study see the following publication:

A.Rafati et al. Chemical and spatial analysis of protein loaded PLGA microspheres for drug delivery applications. Journal of Controlled Release. 162 (2012) 321-329



- o Martyn Davies
- Morgan Alexander
- David Scurr
- o Clive Roberts
- Xinyong Chen



o Alex Shard



- Kevin Shakesheff
- Asme Boussahel



- Shellie Rigby-Singleton
- o Paul Whiteside





- We hope the information provided in this case study is of interest.
- If you wish to get in touch with us to discuss any of the information provided, raise a query/concern or provide feedback then please use any of the methods listed below:

nmRC Commercial Services Nanoscale & Microscale Research Centre University Park Nottingham NG7 2RD

Telephone:+44(0)115 951 5046Email: nmcs@nottingham.ac.ukFax:+44 (0)115 846 7969Website:www.nottingham.ac.uk/nmrc-commercial