



**University of  
Nottingham**

UK | CHINA | MALAYSIA



**University of  
Nottingham**  
Nanoscale and Microscale Research Centre

# **nmRC CASE STUDY**

**EFFECT OF EXCIPIENTS ON SALT  
DISPROPORTIONATION DURING  
DISSOLUTION**

**nmRC\_CS\_03**





University of  
**Nottingham**

UK | CHINA | MALAYSIA



University of  
**Nottingham**  
Nanoscale and Microscale Research Centre

# Effect of Excipients on Salt Disproportionation during Dissolution


## In Situ Raman Imaging Case Study



- An estimated 50-70% of all developed small molecule drugs are administered as salts.
- There is a tendency for the salt to convert back to its free (unionized) form under certain conditions via a reaction known as salt disproportionation.
- The conversion of the salt to the free form is undesirable since it may have detrimental effects on solid-state properties and pharmaceutical product performance.
- In order to understand the dissolution behaviour of pharmaceutical formulations we need to be able study them *in situ* (*i.e.*, during dissolution).




- Oral solid dosage forms are the most common way in which drugs are administered, so dissolution of the dosage form after it is swallowed, specifically the rate at which the active ingredient is released into the body, is of critical importance drug development.
- Standardised dissolution apparatus often rely upon UV-vis absorbance spectroscopy to monitor the drug release profile.






- Raman spectroscopy allows for real-time, '*in-situ*' chemical and morphological analysis that can be used to monitor dissolution dynamics.
- The technique uses the inelastic scattering of light to generate spectra unique to a material's molecular composition and state, and can be used to track changes with time.



***The HORIBA LabRAM HR Raman Microscope available at the Nanoscale and Microscale Research Centre (nmRC), University of Nottingham.***


- In our study, Raman spectroscopy was used to probe the chemical composition of the tablet, whilst UV-vis spectroscopy enabled quantification of the drug release. In order to monitor both simultaneously, a novel flow cell design was used.




*(A) Schematic diagram of the Raman ultraviolet-visible flow cell system. (B) Front view of the flow cell under the Raman spectroscopy objective. (C) Top view of the tablet placed inside the flow cell.*



- The challenge we have applied our novel approach to is the disproportionation of Pioglitazone HCl (PioHCl), which has pH-dependent solubility and therefore release characteristics.




*Pioglitazone HCl*



*pH-solubility profile of a weak base drug presenting the  $pH_{max}$  point and indicating that the solubility may be expressed by two independent regions.*

- Raman spectroscopy enables diagnosis of the composition of the tablet based on the molecular fingerprint of the components.




*Raman spectra of the formulation components within the range 1200–1800  $\text{cm}^{-1}$ .*






# IN SITU ANALYSIS UNDER DIFFERENT CONDITIONS



*Raman spectroscopy: mapping of PioHCl-only tablet during dissolution in (A) acidic and (B) basic media, as a function of time, generated by MCR analysis.*

*UV-vis spectroscopy: dissolution profiles in (A) acidic and (B) basic media, as a function of time, for the PioHCl-only and different formulations.*





- We have demonstrated the utilisation of Raman spectroscopic imaging combined with a UV-vis spectroscopy flow cell is a valuable tool for the investigation of problematic multicomponent pharmaceutical formulations, specifically where disproportionation kinetics depend on the pH microenvironment.
- This approach can enhance our understanding of the interplay between drug and excipient properties, and their influence on disproportionation kinetics, ultimately, helping to mitigate disproportionation of API salts in tablets and enabling pharmaceutical scientists to develop more reliable products with improved drug stability.



**For more details on the work showcased in this case study see the following publications:**

*A. Abouselo et al., "Effect of Excipients on Salt Disproportionation during Dissolution: A Novel Application of In Situ Raman Imaging," Molecular Pharmaceutics, vol. 18, no. 9, pp. 3247-3259, 2021/09/06 2021, doi: 10.1021/acs.molpharmaceut.1c00119.*

The Raman spectroscopy analysis documented here was performed at the Nanoscale and Microscale Research Centre (nmRC) at the University of Nottingham. [www.nottingham.ac.uk/nmrc](http://www.nottingham.ac.uk/nmrc)







- We hope the information provided in this case study is of interest.
- If you wish to get in touch with us to discuss any of the information provided, raise a query/concern or provide feedback then please use any of the methods listed below:

**nmRC Commercial Services**  
**Nanoscale & Microscale Research Centre**  
**University Park**  
**Nottingham**  
**NG7 2RD**

**Telephone:** +44(0)115 951 5046  
**Email:** [nmcs@nottingham.ac.uk](mailto:nmcs@nottingham.ac.uk)  
**Fax:** +44 (0)115 846 7969  
**Website:** [www.nottingham.ac.uk/nmrc-commercial](http://www.nottingham.ac.uk/nmrc-commercial)