



**University of
Nottingham**

UK | CHINA | MALAYSIA



**University of
Nottingham**

Nanoscale and Microscale Research Centre

nmRC CASE STUDY

POLYMER COMPOSITION AND SPATIAL DISTRIBUTION

nmRC_CS_09



University of
Nottingham

UK | CHINA | MALAYSIA



University of
Nottingham
Nanoscale and Microscale Research Centre

Polymer Composition and Spatial Distribution

X-ray Photoelectron Spectroscopy (XPS) Case Study




- Polymers are not always directly or easily identifiable at the surfaces of materials.
- Macro distributions and thicknesses may be of interest, but how can they be assessed?
- X-ray photoelectron spectroscopy (XPS) is a means of chemical analysis sensitive to light elements.
- XPS is also a highly surface sensitive technique (\sim top 10nm), ideal for studying thin films of polymers or liquids*.



* Analyse ionic liquids with the liquid phase photoelectron spectroscopy instrument (LiPPS) at the NNNC.




- X-rays excite electrons out of the sample surface.
- These are collected and their energy analysed.
- Electron energies depend on the element they originated from.
- The resulting spectrum is dependent on what elements are in the surface layer and in what abundance.







- XPS Identification of polymers possible by elemental spectra.
- Example analysis of an ionic liquid:



	Fluorine	Oxygen	Nitrogen	Carbon	Sulphur	Total
No Atoms	6	4	3	14	2	29
At% Expected	20.7	13.8	10.3	48.3	6.9	100



- Parallel imaging of a sample surface can be performed.
- This can provide the spatial distribution of elements, and therefore surface components e.g. polymers
- Example XPS mapping of ionic liquid droplet distribution on a gold surface.





- XPS is a highly surface sensitive technique capable of chemical identification localisation, and quantification.
- Sensitive to light elements with $\sim 0.1\%$ atomic sensitivity it can play a key role in thin film polymer and liquid characterisation.
- Potential applications:
 - ❖ *Polymer identification.*
 - ❖ *Surface enrichments of one polymer constituent.*
 - ❖ *Drug entrapment and depth and coverage with protective layers.*
 - ❖ *Micron scale polymer spatial segregation or coating distributions.*
 - ❖ *Chemical stability of coatings.*
 - ❖ *Polymer mixing in cross section.*



For more details on the work showcased in this case study see the following publications:

E. F. Smith, D. Briggs, and N. Fairley, “Further developments in quantitative X-ray photoelectron spectromicroscopy: preliminary results from the study of germanium corrosion,” *Surf. Interface Anal.*, vol. 38, no. 2, pp. 69–75, 2006.

E. F. Smith, F. J. M. Rutten, I. J. Villar-Garcia, D. Briggs, and P. Licence, “Ionic Liquids in Vacuo: Analysis of Liquid Surfaces Using Ultra-High-Vacuum Techniques,” *Langmuir*, vol. 22, no. 22, pp. 9386–9392, Oct. 2006.

The X-ray photoelectron spectroscopy (XPS) analysis documented here was performed at the Nanoscale and Microscale Research Centre (nmRC) at the University of Nottingham. www.nottingham.ac.uk/nmrc





- We hope the information provided in this case study is of interest.
- If you wish to get in touch with us to discuss any of the information provided, raise a query/concern or provide feedback then please use any of the methods listed below:

nmRC Commercial Services
Nanoscale & Microscale Research Centre
University Park
Nottingham
NG7 2RD

Telephone: +44(0)115 951 5046
Email: nmcs@nottingham.ac.uk
Fax: +44 (0)115 846 7969
Website: www.nottingham.ac.uk/nmrc-commercial