



**University of
Nottingham**

UK | CHINA | MALAYSIA



**University of
Nottingham**

Nanoscale and Microscale Research Centre

nmRC CASE STUDY

**HIGHLY UNIFORM POROUS
MAGNETIC MICROSPHERES FOR
BIOMEDICAL APPLICATIONS**

nmRC_CS_11



University of
Nottingham

UK | CHINA | MALAYSIA




University of
Nottingham
Nanoscale and Microscale Research Centre

Porous Magnetic Microspheres for Biomedical Applications

Mineral Liberation Analysis (MLA) Case Study



- Magnetic nano- and micro-materials have healthcare applications, for example as contrast agents in MRI, magnetically targeted drug-delivery systems, and magnetic-induced hyperthermia treatments.
- Porous microspheres have separate applications in healthcare *e.g.*, drug-delivery systems, tissue regeneration scaffolds etc.
- These materials present functional advantages over dense microspheres *e.g.*, enable payloads to be incorporated on their surfaces or within the pores to enable tailored release kinetics.
- The capacity to develop porous, magnetic materials with high levels of controlled uniformity is highly desirable.





- Porous microspheres are composed of external pores on the surface and internal pores in the core (usually interconnected). Therefore, it is important to explore details of such structures.
- Compositional analysis of such porous micro-systems is crucial to understand **mineral distribution** and **chemical concentration**.
- Important insights into the **chemical reaction pathways** associated with the manufacture of porous microspheres can also be achieved.
- **To investigate homogeneity and internal porosity we need a technique capable of both mineral quantification and high-resolution imaging from a large population of microspheres.**



- The MLA system integrates a combination of back-scattered electron (BSE) imaging and energy dispersive x-ray spectroscopy (EDS) analysis.
- Via software automation this allows visualisation and elemental analysis of a large number of particles for a statistically significant data set.
- Outputs include calculated assay, elemental distribution, particle size distribution, modal mineralogy, mineral association, mineral locking etc.
- Process can take up to 14 polished mounts at a time, or samples up to 150mm² with automated analysis of each sample in turn.
- Resolution reported at 3.5nm @ 30kV





**FEI Quanta 600 with Mineral Liberation Analyser (MLA)
available at the Nanoscale and Microscale Research
Centre (nmRC), University of Nottingham**



Fe_3O_4 /porogen powder

Flame Spheroidisation
process




Sectioning and sieving



Schematic diagram for porous and dense microsphere manufacturing via a rapid flame spheroidisation method, and sample preparation prior MLA analysis






- Secondary electron images reveal **(a)** small, dense and large, porous microspheres; **(b)** localised morphology and topology of porous microspheres; **(c)** microsphere interconnected porosity (boxed in region **(b)**)





(a) BSE image and **(b)** MLA compositional analysis of microspheres, following sieving and sectioning, illustrating microsphere porosity and revealing high levels of homogeneity





Full MLA compositional analysis of FS-processed microspheres, following sieving and sectioning, demonstrating very high levels of $\text{Ca}_2\text{Fe}_2\text{O}_5$

	Name	Pixels	Particles	Area %	Weight %
	Magnetite	318	5	0.01 %	0.01 %
	Magnetite with Ca	138	2	0.01 %	0.01 %
	1 CaFeO	325	3	0.01 %	0.01 %
	2 CaFeO	1866	24	0.07 %	0.08 %
	Ca ₂ Fe ₂ O ₅ (srebr...)	2503262	1501	99.40 %	99.74 %
	4 CaFeO	5481	25	0.22 %	0.14 %
	Calcite	0	0	0.00 %	0.00 %

- A total of **1501 particles (99.74 wt% of the sample)** were classified as $\text{Ca}_2\text{Fe}_2\text{O}_5$ (srebrodolskite).
- EDS data from the MLA showed the elemental composition as follows:

Ca – 19.1%, Fe – 57%, O – 23.9%.



- Mineral Liberation Analysis (MLA) [JKTech/FEI] allows automated large area analysis of polished specimens to identify and quantify mineral composition and distribution.
- Electron microscopy reveals morphologies and porosity of flame-spheroidised porous microspheres.
- MLA technique details:
 - ❖ **Mineral mapping** – used to reveal compositional uniformity of microspheres
 - ❖ **Large number of particles** – used to obtain a statistically significant set of data
- This study demonstrates the role MLA can play in understanding the morphological homogeneity and mineral distribution of microparticles, critical when validating new formulations / manufacturing methods.



For more details on the work showcased in this case study see the following publications:

J. Molinar Díaz, S. A. Samad, E. Steer, N. Neate, H. Constantin, M. T. Islam, P. D. Brown and I. Ahmed,
“Flame spheroidisation of dense and porous $\text{Ca}_2\text{Fe}_2\text{O}_5$ microspheres”, Materials Advances, 2020,
DOI: 10.1039/D0MA00564A





- We hope the information provided in this case study is of interest.
- If you wish to get in touch with us to discuss any of the information provided, raise a query/concern or provide feedback then please use any of the methods listed below:

nmRC Commercial Services
Nanoscale & Microscale Research Centre
University Park
Nottingham
NG7 2RD

Telephone: +44(0)115 951 5046
Email: nmcs@nottingham.ac.uk
Fax: +44 (0)115 846 7969
Website: www.nottingham.ac.uk/nmrc-commercial