

MRC

Medical
Research
Council



MANCHESTER
1824

The University of Manchester

Physically Functionalised Gellan Gum Hydrogels for Neural Cell Culture Applications


Sam Moxon & Nigel Hooper

*Division of Neuroscience and Experimental Psychology,
University of Manchester, UK*

Introduction



- Dementia is now the leading cause of death in England and Wales
- Consequently, research into probing the mechanisms behind dementia using neural cell culture is widespread
- Progress has been made but is often limited by an inability to accurately model the disease *in vitro* as current methods often employ 'classic' 2D monolayer culture on plastic surfaces
- Such strategies are high throughput but don't allow for accurate replication of the soft, complex microenvironment in which neurodegenerative diseases progress
- This can have a profound effect on cell behaviour and inhibit the ability to replicate disease pathologies and interrogate the underlying mechanisms

The Extracellular Matrix



Replicating the ECM *in vitro*

- To overcome issues with monolayer culture, groups often culture cells using *in vitro* 3D cell niches
- Hydrogels are popular tools for this application due to a high water content, low cytotoxicity and ability to replicate microstructural elements of native ECM
- Many food hydrocolloids have been investigated for such applications (gellan gum, alginate, chitosan etc.)



100 μm Stenger *et.al.* (2001)


Collagen-Gellan Blended Hydrogels for Neural Culture

- Brain ECM is mostly comprised of linear, unbranched, gel-forming polysaccharides
- Fibrous proteins penetrate the polysaccharide network allowing for cell-matrix binding
- Blended hydrogels of gellan gum and collagen are being investigated as a simple method to replicate this microenvironment

Collagen-Gellan Blended Hydrogels for Neural Culture

Gellan Gum


- Linear, unbranched polysaccharide of bacterial origin
- Capable of forming porous hydrogels through ionic interactions
- Could replicate the polysaccharide network in brain ECM



Osmalek *et. al.* (2014)



Collagen

- Fibrous, helical protein found in native ECM
- Contains cellular integrin binding domains
- Could replicate the fibrous element of brain ECM responsible for cell attachment



Connery (2015)

Fabrication of Blended Gels




Proof of Concept


- As an initial proof of concept 3 hydrogel blends were prepared with varying gellan:collagen ratios (1:1, 2:1 and 4:1)
- Both polymers were dispersed at a concentration of 0.5% w/w prior to blending
- The final blends were suspended in particulate gel beds of 0.5% agarose
- Frequency sweeps were conducted to determine the mechanical properties of each blend
- SH-SY5Y neuroblastoma cells were also encapsulated in each blend and cultured for 21 days before a live/dead assay was used to determine cell viability and morphology

Rheology

- The storage modulus of brain ECM is reported as being between 0.1 kPa and 1 kPa
- All blends exhibited moduli within this range at low oscillatory frequencies



Live/Dead Assay




Optimising the Microenvironment

- Live/dead images suggested incorporation of collagen improved cell attachment, proving the concept of functionalising gellan with collagen fibres
- However, cells seemed to adhere in clumps suggesting fibres did not penetrate ubiquitously through the gellan network
- Perhaps gelation times were an issue?

Optimising the Microenvironment

- 2:1 and 4:1 ratios of gellan and collagen were prepared suspended in particulate gel beds
- Blends were then incubated at 37°C for 3 different time periods (1 hour, 3 hours and 6 hours) before addition of 100 mM CaCl₂
- This allowed the collagenous component to organise into fibrils for different lengths of time prior to gelation of gellan
- Samples were then freeze-dried and analysed with SEM to analyse the resulting gel networks

SEM Results – 0.5% Gellan



HV	spot	WD	mag	det	Scale
10.00 kV	3.5	13.5 mm	100 x	BSED	500 μ m


0.5% Gellan




HV	spot	WD	mag	det	Scale
10.00 kV	3.5	13.5 mm	500 x	BSED	100 μ m

0.5% Gellan

SEM Results – 0.5% Collagen




HV	spot	WD	mag	det	Scale
10.00 kV	3.5	13.8 mm	100 x	BSED	500 μ m
0.5% collagen					




HV	spot	WD	mag	det	Scale
10.00 kV	3.5	13.6 mm	501 x	BSED	100 μ m
0.5% collagen					

SEM Results – 4:1 (Gellan:Collagen)


1 Hour Incubation




3 Hour Incubation




6 Hour Incubation




1 Hour Incubation



3 Hour Incubation




6 Hour Incubation




SEM Results – 2:1 (Gellan:Collagen)


1 Hour Incubation




3 Hour Incubation




6 Hour Incubation




1 Hour Incubation



3 Hour Incubation



6 Hour Incubation



Current Conclusions/On-Going Work

- Blending 0.5% gellan with 0.5% collagen creates a soft matrix that can better reflect the mechanical properties of brain ECM.
- Live/dead results suggests cells can adhere to blended gels but not to unmodified gellan. Further optimisation is required to achieve more uniform cell adhesion within the whole matrix.
- SEM images could potentially demonstrate increasing incubation time prior to gellan crosslinking results in better formation of collagen fibrils within the matrix but this is difficult to quantify.
- The proof may be in the cell behaviour with the effect of 1, 3 and 6 hour incubation on actin expression now being investigated. This could give us a more quantitative analysis of how cell adhesion can be optimised in the blended gels.
- Once optimised, blended gels will be evaluated for neural cell culture using iPSC-derived neurons.

Acknowledgements

MANCHESTER
1824

The University of Manchester

'Team Hooper'



University of
HUDDERSFIELD

Dr Alan Smith

Nicholas Pearman

Jessica Senior

James Rooney



Manchester Institute
of Biotechnology

MRC

Medical
Research
Council

**Dr Donald
Dean Fund**