



MRI for the Non-Expert

Diffusion MRI

JP Vallée, MD, PhD
Radiology clinics
Diagnostic department
Geneva University Hospital
6

Outlook

Introduction

Principles of diffusion MRI

Data acquisition

Data analysis

Unresolved issues

Outlook

Introduction


Principles of diffusion MRI

Data acquisition


Data analysis

Unresolved issues

DWI earliest marker of stroke



A



B

Fig 2. (A) Twelve axial diffusion-weighted magnetic resonance slices at the maximum b value, showing a large left-sided hyperintense area (volume = 197 cm^3) corresponding to an extensive middle cerebral artery infarction. The National Institutes of Health Stroke Scale Score (NIHSSS) was 20. (B) Axial diffusion-weighted magnetic resonance slice at the maximum b value, showing a small, deep-seated, right-sided hyperintensity (volume = 0.8 cm^3) at the junction of the basal ganglia and the posterior limb of the internal capsule. The NIHSSS was 7.

Prostate cancer




FIGURE 1: A 67-year-old patient, serum PSA 17.8 ng/mL. A low ADC lesion is visible at the mid-level of the prostate. The lesion presents a progressive signal increase on the b-value images starting from low intensity at b0, when compared with the prostatic parenchyma. The low ADC combined with a “flashing” high signal intensity on the b1500 image is typical for a cancerous lesion (arrow). Combining of the ADC map with the b-value images allows for a better accuracy in cancer detection.

Prostate cancer

PI-RADS Assessment


Peripheral Zone (PZ)

DWI	T2W	DCE	PI-RADS
1	Any*	Any	1
2	Any	Any	2
3	Any	–	3
		+	4
4	Any	Any	4
5	Any	Any	5


* "Any" indicates 1-5

What is measured in diffusion MRI (DWI)?


- Displacement of water molecules:
- Typical range 1–17 μm (extracellular)
- Brownian motion, microcirculation....
- Probe of the interstitial space



Reduced diffusion in renal fibrosis




High-grade fibrosis 80%



Normal Absence of fibrosis

$1300 \times 10^{-6} \text{mm}^2/\text{s}$



$2000 \times 10^{-6} \text{mm}^2/\text{s}$

ADC

Outlook

Introduction


Principles of diffusion MRI

Data acquisition

Data analysis

Unresolved issues

Spin motion in a magnetic field




Larmor frequency

Frequency of precession


$$f = \gamma B_0$$

Gyromagnetic ratio

Main magnetic field strength



The diagram shows the equation $f = \gamma B_0$ in the center. Three arrows point from descriptive text to parts of the equation: one from 'Frequency of precession' to the variable f , one from 'Gyromagnetic ratio' to the Greek letter γ , and one from 'Main magnetic field strength' to the variable B_0 .

T2 relaxation





Gradient echo

No B_0 inhomogeneity compensation




B_0 Gradient



B_0 Gradient
inversion






Echoes formation



Modified from Pooley R A Radiographics 2005;25:1087-1099

How to image diffusion?



Outlook

Introduction

Principles of diffusion MRI

Data acquisition

Data analysis

Unresolved issues



DWI acquisition

b0

b20

b300

b700 s/mm²



Multiple acquisitions
Incremental magnetic gradients
(higher gradient = smaller motion probe)
Usually average of 3 spatial directions

DWI acquisition


Table 1. DWI in the kidney: key aspects

Patient preparation	
Hydration	Potential confounder Control by hydrating the patient whenever possible
Data acquisition	
Echo time	Minimize to optimize SNR Minimum limited by maximum b -value
Repetition time	Long enough to allow for T1 relaxation (>1500 ms) Minimum limited by number of image slices
Image orientation	Axial: less motion in image plane Coronal: easier full kidney coverage
Field of view	Usually covers the entire abdomen (320–400 mm) Z-dimension dependent on image orientation
Resolution	Increase: sharpness \uparrow partial volume effects \downarrow Decrease: ETL \downarrow kidney coverage \uparrow SNR \uparrow
ETL	Shorten to lessen susceptibility artefacts Measures: parallel imaging, multi-shot EPI, partial Fourier
Motion compensation	Physiological triggering using external devices Intrinsic triggering using MRI signal (navigator)
b -values	Tailor to respective DWI biomarker Increase number to improve parameter estimates
Image post-processing	
Image quality control	Discard problematic image(s) to ensure imaging parameter value reliability
Motion correction	To account for motion artefacts and eddy current-induced deformations
ROI definition (kidney/medulla/cortex)	From more than one section to have representative average values; no vessels, artefacts, lesions
Model fitting	To compute DWI biomarkers by fitting appropriate signal attenuation models


EPI, echo-planar imaging; ETL, echo train length; MRI, magnetic resonance imaging; SNR, signal to noise ratio

K-space sampling strategy


Conventional K space encoding




Single-shot EPI



Readout-Segmented EPI




(1) Optimization of the acquisition strategy with RESOLVE as a new renal DWI



Friedli I et al MRI 2015

RESOLVE enhanced significantly the quality of renal DWI by reducing susceptibility effects and increasing sharpness.



Outlook

Introduction

Principles of diffusion MRI

Data acquisition

Data analysis

Unresolved issues


ADC (monoexponential fit)

b0



b20

b300


b700 s/mm²




ADC (Coefficient Apparent de Diffusion) [mm²/s]




How to read the images?



Water restriction = DWI \uparrow (high b) = ADC \downarrow




ADC & renal Bx



Berchtold, et al, NDT, 2019

Negative correlation of ADC with interstitial fibrosis in CKD well validated

IVIM fitting



$$S(b) = S_0 e^{(-b \times ADC)} \longrightarrow S(b) = S_b \left((1 - f) e^{(-D \times b)} + f e^{(-(D^* + D) \times b)} \right)$$

Diffusion tensor image (DTI)

Require multiple directions acquisitions (6 - 12)

Parameters:

Fractional anisotropy (FA)

Mean diffusivity (MD)

Tractography




Table 2. DWI biomarker estimation models to investigate renal tissue microstructure

Model	Biomarker(s)	Pros (+) and cons (–)
Monoexponential	ADC: apparent diffusion in the tissue	<ul style="list-style-type: none"> + Most robust against noise + Wide availability and ease of use of biomarker estimation tools – Provides limited information (apparent diffusion only) – Fits DWI data the least
IVIM	D: water diffusion in the tissue D*: pseudodiffusion F: flowing fraction	<ul style="list-style-type: none"> + Describes DWI signal attenuation at best provided sufficient signal-to-noise + Can separate diffusion from pseudodiffusion – No standardized algorithm to compute IVIM parameters
DTI	FA: fractional anisotropy diffusion anisotropy imposed by the tissue microstructure MD: anisotropy-independent mean diffusivity	<ul style="list-style-type: none"> + Provides information on tissue anisotropy – Requires a dedicated acquisition sequence (DTI) with multiple directions
Extended IVIM	D: water diffusion in the tissue D*: pseudodiffusion F: flowing fraction Additional model-specific biomarkers	<ul style="list-style-type: none"> + Potentially advances the characterization of the renal microstructure and microcirculation – Requires complex biomarker estimation – Need further investigation, especially in pathological kidneys
Non-Gaussian	ADC: apparent diffusion in the tissue K/ σ / δ : measure of the degree of deviation of diffusion from a Gaussian law	<ul style="list-style-type: none"> + Accounts for the complexity of diffusion in the renal tissue – Requires complex biomarker estimation – Fits DWI data better than monoexponential but worse than IVIM model

Outlook

Introduction

Principles of diffusion MRI

Data acquisition

Data analysis

Unresolved issues

Unresolved issues

- Best MR sequence / protocol to use (SS-EPI, Resolve, DTI)
- Best analysis to use (ADC, IVIM, DTI)

=> no gold standard

=> Standardisation versus technical evolution!

Unresolved issues

- Origin of the diffusion signal decrease in CKD
(water motion restriction, reduced GFR, reduced perfusion or tubular exchanges...)

=> specificity of DWI ??

but strong validation
(negative correlation of ADC with interstitial fibrosis in CKD)

Unresolved issues

- Validation ++ but qualification is missing
 - Bx replacement in selected CKD
 - New therapy monitoring
 - Selection of Bx patients in acute graft dysfunction
(*Steiger, Eur Rx 2017*)

=> multiparametric & multicentric effort

Acknowledgments

Geneva

Pr S. De Seigneux

Dr L. Berchtold

Dr L. Crowe

Dr I. Friedli

I. Aslam

Fribourg

Pr H. Thoeny

Lausanne

Dr M. Pruijm